






# Lifetime Upgrade Kit for Zeiss LSM-510

- Fluorescence lifetime imaging with ps resolution
- Two-photon excitation
- Non descanned detection
- Ultra-high sensitivity due to TCSPC technique
- Zoom and image rotation functions of LSM-510 applicable
- FRET experiments
- Fluorescence quenching
- New fluorescence probes based on lifetime changes
- Separation of fluorophores by lifetime
- Autofluorescence



The setup employs an advanced multi-dimensional TCSPC technique featuring high count rate, near-ideal counting efficiency, and multi-wavelength capability. It contains the usual building blocks (CFDs, TAC, ADC) in the 'reversed start-stop' configuration together with a scanning interface and a large histogram memory integrated on a single PC board. For each photon the TCSPC module determines the time within the fluorescence decay function and the location within the scanning area. These values are used to address a histogram memory in which the events are accumulated. Thus, in the memory the distribution of the photon density over X, Y, and the time within the fluorescence decay function builds up. The result can be interpreted as a two-dimensional (X, Y) array of fluorescence decay curves or as a sequence of fluorescence images for different times (t) after the excitation pulse.



**Requirements to LSM-510 Microscope**  
 LSM-510 NLO with Ti:Sa Laser  
 Axiocvert or Axioplan version  
 NDD switch box installed




**Becker & Hickl GmbH**  
 Nahmitzer Damm 30  
 12277 Berlin, Germany  
 Tel. +49 30 787 56 32 Fax +49 30 787 57 34  
 info@becker-hickl.com www.becker-hickl.com

US Representative:  
**Boston Electronics Corp**  
 91 Boylston Street, Brookline,  
 Massachusetts 02445 USA  
 Tel: (800) 347 5445 or (617) 566 3821, Fax: (617) 731 0935  
 www.boselec.com tcspc@boselec.com



# LSM 510 NLO FLIM Upgrade Kit



SPC-830 TCSPC module



DCC-100 detector controller

## Detector Options



### PMC-100 / M-SHUT-Z-NDD

Detector / shutter assembly  
Cooled PMT  
IRF width 150 to 180 ps  
For lifetimes down to 100 ps  
Dark counts typ. 20 to 50 s<sup>-1</sup>  
PMC-0: 300 to 650 nm  
PMC-1: 300 to 820 nm  
Length 135 mm<sup>1)</sup>  
Gain control via DCC-100  
Overload shutdown via DCC-100  
no external components



### H7422-40 / M-SHUT-Z-NDD

Detector / shutter assembly  
Cooled GaAsP PMT  
IRF width 300 to 350 ps  
For lifetimes down to 200 ps  
Dark counts typ. 20 to 80 s<sup>-1</sup>  
300 to 690 nm, ultra-high sensitivity

Length without cables 195 mm<sup>1)</sup>  
Gain control via DCC-100  
Overload shutdown via DCC-100  
HFAC-26-2 preamplifier



### R3809U / M-SHUT-Z-NDD

Detector / shutter assembly  
Multichannel plate PMT  
IRF width < 30 ps  
For lifetimes down to 20 ps  
dark counts typ. 100 to 500 s<sup>-1</sup>  
R3809U-52: 180 to 650 nm  
R3809U-50: 180 to 820 nm  
Length without cables 195 mm<sup>1)</sup>  
Gain control via DCC-100  
Overload shutdown via DCC-100  
HFAD-26-01 preamplifier  
FuG HCN-14-3500A HV power supply

1) Depending on the microscope configuration, space constraints may preclude the use at one of the detector ports of the NDD switch box

## Dual Wavelength Operation




### HRT-41 Router

Simultaneous operation of two PMC-100, H7422, or R3809U detectors attached to the outputs of the Zeiss NDD switch box

Uses BH's patented simultaneously recording multi-detector TCSPC principle

## FLIM Data Analysis



### SPCImage Data Analysis

Single and multi exponential decay analysis  
Lifetime in pixels displayed as colour  
Lifetime distribution  
FRET intensity by double-exponential analysis of donor decay function

## For more information please download or call for

SPC-134 through SPC-830 TCSPC modules, Manual and TCSPC compendium  
Upgrading TCSPC laser scanning microscopes with the SPC-730 and SPC-830 TCSPC imaging modules  
W. Becker, A. Bergmann, Lifetime imaging techniques for optical microscopy  
W. Becker, A. Bergmann, Detector assemblies for picosecond microscopy  
DCC-100 manual  
SPCImage manual  
SPC-830, DCC-100, PMC-100, HFAC-26 and HFAD-26 data sheets

### Becker & Hickl GmbH

Tel. +49 30 787 56 32 Fax +49 30 787 57 34  
info@becker-hickl.com www.becker-hickl.com

