
PML-Spec

Multi-Spectral Lifetime Detection


Multi-wavelength detection of fluorescence decay functions


16 simultaneously recording wavelength channels
Spectral range 300-850 nm
High time resolution: 180 ps fwhm IRF width
Useful count rate > 2 MHz
Ultra-high sensitivity
Short acquisition times
Greatly reduced pile-up
Works with all bh TCSPC modules


Biomedical fluorescence
Autofluorescence of tissue
Time-resolved laser scanning microscopy
Multi-spectral lifetime imaging
Recording of chlorophyll transients
Stopped flow fluorescence experiments



The PML-SPEC uses bh's proprietary multi-dimensional TCSPC technique. The light is split into its spectrum by a polychromator. The spectrum is detected by a 16-channel multi-anode PMT. The single photons detected in the PMT channels are recorded in a bh TCSPC module. The TCSPC module builds up a photon distribution over the time in the fluorescence decay and the wavelength. The technique does not use any time gating, detector channel multiplexing, or wavelength scanning and therefore reaches a near-ideal counting efficiency.










Becker & Hickl GmbH Nahmitzer Damm 30 12277 Berlin, Berlin Tel. +49 / 30 / 787 56 32 Fax. +49 / 30 / 787 57 34 email: info@becker-hickl.com www.becker-hickl.com



US Representative: Boston Electronics Corp tcspc@boselec.com www.boselec.com



UK Representative: Photonic Solutions PLC sales@psplc.com www.psplc.com



Covered by patent DE 43 39 787

PML-Spec

Multi-Spectral Lifetime Detection

Specification

Optical System

Input slit height, mm

Type of grating, lines/mm	400	600	1200
Recorded interval ¹ , nm	320	208	106
Wavelength channel width, nm	20	13	6.65
Spectral range of grating ² , nm	300-850	300-850	300-850
F number		F/3.7	
Input slit width, mm		0.6	


¹ any interval within spectral range of grating

Detector³

Detector		
Cathode spectral response	bi-alkali, 300 to 600 nm	multi-alkali, 300 to 850 nm
Typical dark count rate, s ⁻¹	200	800
Number of spectral channels	16	
Timing output polarity of detector	negative	
Average timing pulse amplitude	40 mV	
Time resolution (FWHM)	150 to 200 ps	
Time skew between channels	< 40 ps	
Timing output connector	SMA, 50Ω	
Routing signal	4 bit + Count Disable Signal, TTL/CMOS	
Routing signal connector	15 pin Sub-D / HD	
Power supply	-800900V / 0.35 mA from external HV power supply	
* * *		1 11 2

3

³ please see data sheet and manual of PML-16 multichannel PMT head



Application

Wavelength-Resolved Fluorescence Decay Measurement

Related Products and Accessories

SPC-134 through SPC-830 TCSPC boards, ps diode lasers, FLIM upgrade kits for scanning microscopes, fibre couplers, computer-controlled HV power suplies. Please see www.becker-hickl.com or call for individual data sheets.



Becker & Hickl GmbH
Nahmitzer Damm 30
12277 Berlin, Berlin
Tel. +49 / 30 / 787 56 32 Fax +49 / 30 / 787 57 34
iwww.becker-hickl.com info@becker-hickl.com

Boston Electronics Corporation

91 Boylston Street, Brookline.

Massachusetts 02445 USA
Tel: (800) 347 5445 or (617) 566 3821, Fax: (617) 731 0935
www.boselec.com tcspc@boselec.com



² limited by spectral response of detector